

MOUNTAIN HIGH OXYGEN

CMK CYLINDER BAND & CUSHION MATRIX

CYLINDER ID	KIT PART NO.	DESCRIPTION	BAND DIAM.	BAND PART NUMBER	KNOB	FOOT	CUSHION
							Flat Length
AL-113/180	00CMK-0010-02	QR	3.25	00CMK-0040-02	Short	Short	NA
AL-248/415/682	00CMK-0016-02	QR	4.50	00CMK-0048-02	Short	Short	NA
AL-647	00CMK-0004-02	QR	5.25	00CMK-0052-02	Long	Short	NA
CFF-480	00CMK-0022-03	QR W/CUSHIONS	4.75	00CMK-0050-02	Short	Short	6"
CFFC-048/KF-50	00CMK-0006-03	QR W/CUSHIONS	7.12	00CMK-0055-02	Long	Tall	9"
	00CMK-0006-13	QR W/CUSHIONS	7.12	00CMK-0055-02	Long	Saddle	9"
CFFC-022	00CMK-0004-03	QR W/CUSHIONS	5.25	00CMK-0052-02	Long	Short	6"
CFFC-025	00CMK-0004-03	(DD-Lite)	5.25	00CMK-0052-02	Long	Short	6"
KF-011	00CMK-0002-03	QR W/CUSHIONS	3.75	00CMK-0044-02	Short	Short	4"
KF-077	00CMK-0008-03	QR W/CUSHIONS	7.75	00CMK-0070-02	Long	Tall	9"
KF-115	00CMK-0012-03	QR W/CUSHIONS	9.25	00CMK-0072-02	Long	Tall	12.5"
	00CMK-0012-13	QR W/CUSHIONS	9.25	00CMK-0072-02	Long	Saddle	12.5"
Foot, Short	00CMK-0030-00						
Foot, Long	00CMK-0032-00						